

Digital Healthcare Innovator

# INTIN

Innovation

Interesting

Information

## INTIN is with you in your everyday life!



**Headquarters / R&D / Manufacturing  
Innovation Center**

South Korea, Daegu Innovation City  
South Korea, Gyeonggi-do, Pangyo

**Phone** +82 53-963-7575

**Email** [intin@intin.kr](mailto:intin@intin.kr)

**Website** [www.intin.kr](http://www.intin.kr)

Are you still hesitant to go to the hospital?

# If you add **INTIN** to your troubles —



Did you know?  
The fact that our bodies sends us a lot of signals before our health condition is lowered or is in bad condition?

INTIN's OVIEW Digital Medical devices and Personal Health Service APPs help solve these kinds of problems and concerns

## Your healthy live begins!

We know it's important, but it's also hard to practice good health.

Digital Healthcare that helps you live a healthy life - OVIEW



Healthy Life Design starts with

## Digital healthcare – OVIEW

Healthy preparation for plans to have children

Safely managed sex life

Preparing for menopause in advance

Analysis of air quality to manage the air quality of the house

Management of respiratory diseases for the whole family

Checking on male health in the comfort of a home

# OVIEW Sperm Tester & OVIEW Sperm Tester Pro

OVIEW Sperm Tester and Sperm tester Pro allows you to check sperm population and activity though self-test at home - solving your concerns about the traditional uncomfortable and embarrassing sperm testing methods.



Oview sperm tester Pro



Oview sperm tester

### | Countries where available



## CLINICAL TRIALS


- Manpower Urology(2018): Clinical trial comparison between Oview Sperm Tester Pro and a semen analysis microscope(200 subjects)
- Pusan University Hospital (2021): Clinical trial comparison between Oview Sperm Tester Pro and a semen analysis microscope(35 subjects, in progress)
- Hanyang University Hospital (2021): Clinical trial comparison between Oview Sperm Tester Pro and CASA(35 subjects, in progress)

99.1% Accuracy

## Manage your health more accurately and conveniently with your smartphone APP!

Record my health with your smartphone. Get more accurate information with AI and Big Data technology.

- ✔ Film and check the results of the OVIEW through your smartphone
- ✔ Advices from Urology, Psychiatry, and Food and Nutruition for sperm quality improvement
- ✔ 12 week health improvement program
- ✔ History management capabilities



\* Search for "OVIEW M" in the Google Playstore / App Store

# OVIEW Sperm Tester Pro

## User Guide

Check out the user guides in the OVIEW Youtube channel.

OVIEW INTIN ▶  
Youtube QR code



**01 Semen Liquefaction**

15 Min.

**02 Insert Semen**

Collect here

**03 Insert Chamber**

**04 Attach to Smartphone**

**05 Start Analysis**

**06 Check Results**

Sperm Count	
High	802 M
Medium	255 M
Low	39 M

Unit: Million/ml

Sperm Motility	
High	78 %
Medium	61 %
Low	40 %

Sperm Count	
High	802 M
Medium	255 M
Low	39 M

Unit: Million/ml

Sperm Motility	
High	78 %
Medium	61 %
Low	40 %

When you want to check the exact ovulation date / menstrual cycle

# Oview Ovulation Tester

OVIEW Ovulation Tester allows you to easily measure your ovulation with saliva anywhere, anytime. Predict accurate ovulation days and menstrual cycles through the smartphone application!



- ✓ Film and check the results of the OVIEW through a Smartphone
- ✓ Oview-Calendar function (Ovulation date / Expected Menstrual Cycle check)
- ✓ Ovulation · menstrual day notification
- ✓ Diary Password
- ✓ Result sharing

| Countries where available



South Korea



Japan



Thailand



Russia

## CLINICAL TRIALS

- Yeungnam University Hospital (2017):
- Hanyang University Hospital, Incheon St.Mary's Hospital, Korea University Guro Hospital (2021):


Clinical trial comparison between Oview Ovulation Tester and urine ovulation analyzer (100 subjects)

95.1% Accuracy

Clinical trial comparison between Oview Ovulation Tester, ultrasound, and blood test. (35 subjects, in progress)

### Fertility Window

- ESTRADIOL
- FSH
- PROGESTERONE
- LH



※ In the ovulation period, when estrogen increases, the salt concentration of saliva increases together, and the crystal structure of saliva changes to the shape of bracken leaves.

# OVIEW Ovulation Tester

## User guide

Check out the user guides in the OVIEW Youtube channel.

OVIEW INTIN ▶  
Youtube QR code



**01 Collect Saliva**

**02 Take care of strip orientation**

**03 Insert strip into the device**

**04 Attach clip with the device**

**05 Analyze**

**06 Check Results**


The diagram illustrates the final step of the process. On the left, a circular field contains several small dots. An arrow points to the middle, where the field is filled with a complex branching structure. A second arrow points to the right, where a smartphone screen displays the results. The screen shows a circular field with the branching structure, and below it, a results summary:

Result	Ovulation
Fertility Status	High
Information	
Pain	Min Middle Max

When you need, non-face-to-face, personal respiratory device

# Oview Multi - Respiratory Diagnosis / treatment device

From diagnosis to treatment of respiratory diseases, the OVIEW Multi - Respiratory diagnosis / treatment device can be used conveniently using various replaceable modules on a single device.



Oview Multi - Respiratory Diagnosis / treatment device

## Use it like this!



Record treatment module usage and manage preventive care through follow-up management in OVIEW Multi APP



Oview indoor Air monitor checks/manages the air quality of your house



# ONESTOP SYSTEM - MULTI

OVIEW INTIN ▶  
Youtube QR code



## Why do you need Oview Multi?



Air Pollution



Smoking



Infection



Allergy

Total respiratory health management

Convenient



Modular



**Regular hospital visits and 'Care' of your health for comfortable breathing in daily life are essential.**

## Air monitor

Check your air environment around you with an indoor air analyzer! After, remove foreign substances from the nose to treat your respiratory system with OVIEW Multi  
(You can choose different products depending on the disease)



# Oview Multi / Oview




## Thermometer

Accurately measure body temperature without contact using infrared sensors.



## Multi Nose Spray

Uses saline, etc. to clean and remove foreign substances from your nose.




## Multi Nebulizer

Converts prescribed medicine into a microscopic gaseous particle for patients with bronchial asthma or chronic fractured lung disease to easily inhale.



## Multi App

Can be linked to the OVIEW Air Monitor to check the air quality and diagnose and treat user's diseases. Easily check and track data that can help you keep track of your health in your daily life.






## Multi Spirometry

Measure your current health by checking your breathing volume.




## Multi Nose Cleaner

Removes foreign substances and mucus from the nasal cavity through suction.




## Multi LED Care

Increases blood flow in the nose to relieve rhinitis using a LED light source.



## OVIEW Air monitor

Information, such as temperature, humidity, fine dust concentration and VOC, are shared and managed through application.



Room temperature, humidity, fine dust density and VOC are measured.

Micro-dust ( $\mu\text{g}/\text{m}^3$ )

Concentration	0-15	16-50	51-100	101-
Icon				
LED Color	<span style="color: blue;">—</span>	<span style="color: green;">—</span>	<span style="color: orange;">—</span>	<span style="color: red;">—</span>

Follow INTIN on Social Media  
for various news about Oview!

Instagram @intin\_oview @intin\_oview\_korea  
Facebook @INTIN OVIEW @INTIN OVIEW 오뷰

Design

your

health

3.5M

o:view  
SMART OVULATION TESTER

INTIN  
INTIN  
INTIN

INTIN  
Health Platform

INTIN Inc. CS CENTER +82-53-963-7575 E intin@intin.kr W intin.kr/en  
A 52, Cheombok-ro, Dong-gu, Daegu, Republic of Korea © 2021 INTIN. All Rights Reserved.